

PHYSIOMER®

All the goodness of the sea to breathe better.

PHYSIOMER® is a complete range of natural seawater nasal care products.

PHYSIOMER® Kids

Nasal Spray

Hygiene - Active prevention

From 2 years old

115 ml

Properties

- Eliminates excess mucus due to cold
- Cleans and moistens nasal passages
- Used regularly, prevents the recurrence of cold and the risk of secondary infections and otitis

100% natural

- Isotonic seawater
- Preservative free
- Suitable for daily use

TO KNOW BEFORE USING PHYSIOMER® Kids

- Use by more than one person may spread infection.
- In case of nasal irritations, temporary tingling may appear during the use.
- Do not use after the expiration date printed on package.
- Do not use if the packaging has lost its original appearance.
- Do not use a broken or damaged product.
- If nasal symptoms prevail beyond 10 days for adults and children, consult a healthcare professional.
- Keep out of reach and sight of children.

INGREDIENTS

100% natural, undiluted and isotonic seawater solution (equivalent to 9 g/l of NaCl). Sterile and preservative free.

HOW TO USE

Recommended use

- In treatment of nasal symptoms, in case of cold (runny and/or stuffy nose) or flu: 4 to 6 times a day, spray 1 to 2 seconds in each nostril.
- In daily use, in hygiene and prevention: 1 to 2 times a day, spray 1 to 2 seconds in each nostril.

Administration mode

Adults and children: Adult supervision is recommended for use in children.

Before use: check if the spray works properly.

- 1 - Tilt the head to one side. Gently insert the nozzle in the upper nostril and press the nozzle once firmly for 1-2 seconds. Let the solution act for a few seconds.
- 2 - Repeat with the other nostril.
- 3 - Blow your nose.



Wash the nozzle with hot water after each use and dry properly.