PHYSIOMER® Baby

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Babies are particularly affected when nasal blockage or excess mucus: they have not yet learnt to blow their nose. Therefore, babies may be more susceptible to complications resulting from a blocked nose, including inflammation of the ears and air passages.

PHYSIOMER® Unidoses gently cleans your baby's nasal passages and decreases susceptibility to ear and respiratory passages infection.

Unidoses

From birth

PROPERTIES:

- Gently cleanses and moisturizes baby's nasal passages.
- Thins and helps remove excess mucus.
- Helps to ease breathing.
- Helps to clear dust, allergens, pollutants and pathogens (virus and bacteria).
- Contributes to prevent the occurrence of ENT complications (such as otitis).

INGREDIENTS:

- 100% natural, isotonic seawater solution (equivalent to 9 g/l of NaCl).
- Sterile and preservative free.

Other products available to relieve your child's nasal passages:

PHYSIOMER® Baby

Nasal Aspirator

From birth

Gentle, effective and quick way to relieve nasal congestion. Equipped with a soft nozzle.

PHYSIOMER® Baby

115 ml

From 15 days of age

Cleans your toddler's nasal passages and helps restore breathing as well as comfort and wellbeing during sleeping and feeding.

PHYSIOMER® Kids

115 ml

From 2 years old

Restores your child's nasal breathing from a stuffy nose and prevents the frequent cold complication into otitis. Daily, it helps protect your child from repetitive nasal infections upon entering daycare or school.

Also convenient for adults for hygiene ocular use. Product availability changes according to countries.

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Break off one dose and turn the tip to open.

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FOR NASAL WASH

For babies and newborns: It is recommended to proceed with nasal wash before feeding, when the child is awake and calm. Lay the baby down on a slightly inclined bed (30°) with his/her head to one side (never backwards to avoid the risk of wrong way into the throat).

1/ Gently insert the tip in the upper nostril and press the single dose firmly to expel the solution within the nostril. If needed, press a second time in order to fully empty the single dose.

2/ Let the solution flow out through the other nostril.

3/ Repeat with the other nostril.

For babies and newborns: lift the baby's back up, and wipe his/her nose.

For better results in case of runny and/or stuffy nose, use preferably with **Physiomer® Nasal Aspirator**.

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FOR EYELID CARE

1/ Wet a sterile compress with PHYSIOMER® Unidoses.

2/ Wipe the eyelids with the compress, make sure to avoid spreading dust, dirt or infection.

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RECOMMENDED POSOLOGY

Nasal use: Can be used as often as needed.

A dose per nostril up to 4 times a day.

Ocular use: Can be used as often as needed.

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Do not re-use an opened dose to prevent any risk of infection or secondary infection. **Not for injection.**

If nasal symptoms prevail beyond 5 days for babies and 10 days for adults and children, consult a healthcare professional.

Use by more than one person may spread infection.

Do not use if the packaging has lost its original appearance.

Do not use after the expiry date printed on package.

Keep out of reach and sight of children.

For babies: do not perform nasal wash after feeding.

Storage conditions should not exceed 27°C.

Predictable adverse effects: In case of nasal irritations, temporary tingling may appear during the use.

Contraindications: No contraindications have been identified.